# **PRACTICAL CLINICAL COURSES**

A Service of the Gordon J. Christensen Career Development Program

# X1946 Conservative Indirect Restorations

Jose-Luis Ruiz, DDS Gordon J. Christensen, DDS, MSD, PhD

**Materials Included:** 

C.E. Instruction Sheet AGD Post-Test

# Gordon J. Christensen PRACTICAL CLINICAL COURSES

# PROCEDURE FOR RECEIVING CE CREDIT FOR VIDEO COURSES

- Complete the enclosed Post-Test.\* For each <u>CE Video Purchased</u>, one complimentary test is included.
   Additional tests are \$35/person per test. Fees can be paid by check or with a credit card when submitting tests to Practical Clinical Courses. If more than 10 CE tests are required per video, please contact us at 800-223-6569.
- 2. Complete the demographic information located at the end of the test.

#### Type of Credit:

- a. If the applicant selects "State," PCC will issue a certificate of verification to the applicant. The applicant must then submit this certificate to their state board to obtain credit.
- b. If the applicant selects "AGD," PCC will submit credit information to the Academy of General Dentistry and confirmation to the applicant that it has been submitted (the applicant may check their AGD transcript for verification please allow 30 days).
- c. If the applicant selects "Both," PCC will complete a. & b. above.
- 3. Return the **Post-Test portion** via mail, fax, or email. Our contact information is as follows:

Practical Clinical Courses 3707 N Canyon Road Suite 3D Provo, UT 84604 Fax: (801) 226-8637 cetests@pccdental.com

4. Practical Clinical Courses will correct the Post-Test. Passing scores are 70% or higher.

\*TO OBTAIN CE CREDIT ONLINE: Login or create an account on <a href="www.pccdental.com">www.pccdental.com</a> and select "My CE Tests" from the left-side menu. Click on the video title to take the test online. RESULTS ARE IMMEDIATE. Missing the test? Contact us at 800-223-6569 during our business hours of 7:00 a.m. – 4:00 p.m. MST to add it to your account.

## **POST-TEST**

#### **X1946 Conservative Indirect Restorations**

1. Why use minimally invasive dentistry?

	<ul><li>a. Increases longevity of natural dentition</li><li>b. Preserves enamel</li><li>c. Is usually less expensive for patients</li><li>d. All the above</li></ul>
2.	Placing an onlay instead of a crown helps preserve the pulp and prevent future root canal therapy.  a. True  b. False
3.	Subgingival crown margins can experience:  a. caries from open margins on crowns.  b. chronic inflammation.  c. receding gingiva.  d. all the above.
4.	Which is the longest-lasting retention feature?  a. Enamel impregnated with resin cement  b. Bond to dentin
5.	<ul> <li>When are onlays a potential alternative?</li> <li>a. When one cusp is broken off</li> <li>b. When a tooth needs restoration after endodontic therapy</li> <li>c. Almost any time a crown is indicated</li> <li>d. All the above</li> </ul>
6.	Dentists can use the enamel margin preservation approach for:  a. direct restorations.  b. indirect restorations.  c. both a. and b.  d. none of the above.
7.	Margin elevation is an alternative for deep-crown margins.  a. True  b. False
8.	On the external surface of a tooth, there are organisms.  a. 150 b. thousands of c. millions of d. billions of

## **POST-TEST**

### **X1946 Conservative Indirect Restorations**

9.	When applying two 1-minute applications of GLUMA or MicroPrime, organisms on the tooth surface will:			
	a.	be killed on the enamel only.		
	b.	be killed on the dentin only.		
	c.	be killed to the pulp.		
	d.	not be killed completely.		
10.	Scanning is more accurate than using vinyl or polyether.			
	a.			
	b.	False		
11.	IPS e.max crowns have three different opacities. Dr. Ruiz suggested using for anterior teeth and for posterior teeth.			
	a.	IPS e.max MT, IPS e.max LT		
	b.	IPS e.max LT, IPS e.max MT		
	c.	IPS e.max MT, IPS e.max HT		
	d.	IPS e.max HT, IPS e.max MT		
12.	. Excellent margin blends need cement with translucency.			
	a.	light		
	b.	medium		
	C.	dark		
	d.	any type of		
13.	When preparing teeth for porcelain veneers, a incisal edge thickness is usually sufficient.			
	a.	0.5 mm		
	b.	1.0 mm		
	c.	1.5 mm		
	d.	2.0 mm		
14.	E.max is one color (one opacity), so variations in thickness of the tooth preparation will not affect the color.			
	a.	True		
	b.	False		
15.	Having a system in place for smile design, such as the dento-facial esthetic diagnosis system, allows you and			
	your pa	itients to assess patient expectations of:		
	a.	color.		
	b.	shape.		
	C.	length.		
	d.	all the above.		
16.	Minimally invasive dentistry is more difficult treatment than conventional dentistry.			
	a.	True		
	b.	False		

#### **POST-TEST**

#### **X1946 Conservative Indirect Restorations**

- 17. Canine guidance is crucial because it decreases occlusal forces.
  - a. True
  - b. False
- 18. Preserving at least 50% of the volume of surface enamel is necessary for adequate veneer retention.
  - a. True
  - b. False
- 19. Some of the benefits for using supra-gingival minimally invasive dentistry include:
  - a. easier impressions, less post-op pain, longer lasting.
  - b. fewer problems if temps fall off, healthier, less chance of RCTs.
  - c. both a. and b.
  - d. none of the above.
- 20. If you place veneers, use a dual-cure cement. If you place onlays, use a light-cure cement.
  - a. True
  - b. False

#### **PLEASE PRINT**

Name				
Name of video purchaser (if different from above)				
Address				
City/State/Zip				
Phone No.				
Email				
Indicate which type of credit you wish to obtain: □ AGD □ State □ Both				
State License No	AGD No			
Payment information   Visa   American Express   Mastercard   Discover				
Card #	Expires/ CVV2-Code:			

For each <u>CE Video Purchased</u>, one complimentary test is included. If you require CE tests for *staff members or additional doctors to receive credit, the fee is \$35/person per test*.

If more than 10 CE tests are required per video, please contact PCC at 800-223-6569.